Supplemental Document

## **Optics Letters**

## Sub-femto-Newton sensing torsion pendulum for detection of light force: supplement

SHENGGUO GUAN,\* JING SUN, BIAO HUANG, YANBEI CHENG, ZHENGLU DUAN, AND JIANXIN LE


College of Physics and Communication Electronics, Jiangxi Normal University, Nanchang 330022, China \*Corresponding author: gsg@jxnu.edu.cn

This supplement published with Optica Publishing Group on 21 September 2022 by The Authors under the terms of the Creative Commons Attribution 4.0 License in the format provided by the authors and unedited. Further distribution of this work must maintain attribution to the author(s) and the published article's title, journal citation, and DOI.

Supplement DOI: https://doi.org/10.6084/m9.figshare.20980459


Parent Article DOI: https://doi.org/10.1364/OL.472130

## Sub-femto-Newton sensing torsion pendulum for detecting of light force: supplemental document




The picture of the experimental device and the key parameters of it are as follows:

Fig. S1. Panoramic picture of experimental device (PC not displayed). The experimental system is placed in the laboratory in the basement. The laboratory has a good vibration isolation and constant temperature environment ( $\Delta T \sim 0.05 \text{ K/day}$ ).



 $\label{eq:Fig.S2} Fig. S2. \quad The signal laser beam power changes by 1.07 mW, and optical attenuators attenuates it to about 0.7 \mu W.$ 



 $Fig.~S3. \quad Angle ~detection ~system. ~The ~position ~detector ~is ~about ~1.31m ~away ~from ~the ~torsion ~pendulum.$ 





Fig. S4. Inside of vacuum chamber. The torsion pendulum is suspend by a 46 cm long tungsten fiber,  $8\mu$ m in diameter, and an aluminum-coated silicon wafer, 70 mm×5 mm×0.2mm in size, which has a mass of 0.17 g, the total mass of torsion pendulum is about 1.4 g.



The torsion pendulum

Fig. S5. Another view of the torsion pendulum. A pair of capacitor plates with an area of 0.6cm<sup>2</sup> are used for dynamic amplitude control, they are placed 2cm opposite the torsion pendulum and 2.5cm away from the tungsten wire.