Optics EXPRESS

Liquid-crystal based drift-free polarization modulators: Part II. Ultra-stable Stokes and Mueller polarimeters: supplement

JEAN REHBINDER,^{*} JEAN DELLINGER, BRISÉIS VARIN, MARC TORZYNSKI, YOSHITATE TAKAKURA, CHRISTIAN HEINRICH, AND JIHAD ZALLAT

Laboratoire ICube, Université de Strasbourg, Bd Sébastien Brant, 67412 Illkirch, France *rehbinder@unistra.fr

This supplement published with Optica Publishing Group on 10 March 2023 by The Authors under the terms of the Creative Commons Attribution 4.0 License in the format provided by the authors and unedited. Further distribution of this work must maintain attribution to the author(s) and the published article's title, journal citation, and DOI.

Supplement DOI: https://doi.org/10.6084/m9.figshare.22110497

Parent Article DOI: https://doi.org/10.1364/OE.480774

Liquid-crystal based drift-free polarization modulators - Part II. Ultra-stable Stokes and Mueller polarimeters: supplemental document

Simulating LCVR operation with an LCM. Figure S1 compares the evolution of the retardance of an LCM in its normal configuration (compensated temperature dependence) and in a setting simulating an LCVR ($\varphi_B=0$).

In this experiment, the temperature in the enclosure was set successively to 22, 38, and 50°C. The retardance of the LCM was measured for a set of ϕ_A and ϕ_B values.

The figure illustrates how the temperature-dependence vanishes at low retardance in the LCVR setting. It also shows the effectiveness of the LCM for obtaining stable retardance values over a wide temperature range.




Fig. S1. Temperature dependence of the retardance of a LCM with compensation (LCM, blue lines) and without compensation (simulating a LCVR, red dotted lines).

Repeatability of the measurements. We compare two measurements with the same system (see section 3 and figure 3) and the same configuration as the experiments described in the main document. Measure 1 and measure 2 were performed at few days' interval.

Figure S2-S4 are counterparts to Figure 4-6 of the main document.




Fig. S2. Temperature dependence of the measurement matrix *A* of a PSA with compensation (LCM-PSA, blue lines) and without compensation (LCVR-PSA, red lines). Full lines of light color correspond to measure 1, dotted lines of darker color to measure 2.




Fig. S3. Temperature dependence of the transfer matrix $A_{25^{\circ}C}$ ⁻¹ A_{Tmeas} of a PSA with compensation (LCM-PSA, blue lines) and without compensation (LCVR-PSA, red dotted lines). Full lines of light color correspond to measure 1, dotted lines of darker color to measure 2.




Fig. S4. Temperature dependent parameters of a waveplate obtained with compensated (LCM-PSA) and standard (LCVR-PSA) PSA in the Mueller configuration: (a) Retardance, (b) orientation of the retarder, (c) Diattenuation, and (d) Depolarization extracted from the Mueller matrix using polar decomposition